

EX00107 Sample preparation protocol

Samples were stored in culture plates (9cm). Per protocol, plate was washed with ice cold MeOH/H₂O/CHCl₃ (80:10:10 by Vol.). Plate was scraped with disposable spatula and rocked to and fro for 20-30 seconds with crash solvent. Solvent removed in it's entirety to a 2mL deep well. Sample order noted for LIMs upload sheet. Deep well was then capped and placed in centrifuge set to 4°C, 4750RPM (5251g) for 60 min. Post spin 400µL of all samples were transferred to a Agilent vial insert for dry down under 100% nitrogen. Once dried all samples but blank (100% H₂O) was reconstituted with 200µL of TBoc (6.25µL of stock standard solution in 1mL of H₂O).