

## LC - Instrument Parameters AbsoluteIDQ® p180 Kit (Biocrates Life Science AG, Innsbruck, Austria)

### System:

UPLC (Waters Acquity, Waters Corporation, Milford, USA)

### Column:

Agilent, Zorbax Eclipse XDB C18, 3.0 x 100 mm, 3.5 µM, Agilent Waldbronn, Germany

### Precolumn:

Security Guard, Phenomenex, C18, 4 x 3 mm; Phenomenex, Aschaffenburg, Germany

### Materials:

Water MiliQ,

Methanol (Merck KGaA, Darmstadt, Germany, hypergrade for LC-MS)

Acetonitril (Merck KGaA, Darmstadt, Germany hypergrade for LC-MS)

Ammonium Acetate (Honeywell - Fluka, Seelze, Germany)

Formic Acid (Honeywell, Fluka, Seelze, Germany)

Running Solvent: 5mM ammonium acetat in methanol

FIA: B: Mix running solvent and MiliQ water 1:1

### Gradient Table FIA

Time (min)	Flow Rate (mL/min)	%A	%B	Curve
Initial	0.030	0.0	100.0	Initial
1.60	0.030	0.0	100.0	6
2.40	0.200	0.0	100.0	6
2.80	0.200	0.0	100.0	6
3.00	0.030	0.0	100.0	6