

## Sample preparation for MetIDQ p180 Kit measurement

### Solvents:

Acetonitril (Merck KGaA, Darmstadt, Germany hypergrade for LC-MS)

Water MiliQ,

Extracting agent - ACN / H<sub>2</sub>O (1:1)

### Equipment

4 steel balls size M

Eppendorf Tubes 2mL

Tissue slicer (Rettberg, Germany)

Centrifuge (Sigma)

### Work steps

Approximately 100 mg of sample and 4 steel balls of size M into each tube. Per mg of sample 5 µL of extracting agent was added. Shake the samples for 10 minutes at 30 Hz in the tissue slicer and centrifuge for 2 minutes at 14000 rpm. 10µL of supernatant was used for the targeted analytics.

For blood sample analysis 10 µL plasma was taken.

### Kit reparation

The analysis was performed using the validated MetIDQ p180 Kit and described in Siskos et al. [1]. Data processing was carried out with the provided quantitation method Kit (Biocrates Life Sciences AG, Innsbruck, Austria).

1 Siskos AP, Jain P, Römisch-Margl W, Bennett M, Achaintre D, Asad Y, *et al.* Interlaboratory Reproducibility of a Targeted Metabolomics Platform for Analysis of Human Serum and Plasma. *Anal Chem* 2017;**89**:656-65.