

GC-MS analysis of tricarboxylic acid cycle (TCA) metabolites

Targeted profiling of TCA metabolites was conducted at the Mayo Clinic by gas chromatography mass spectrometry (GC-MS) as previously described^{55,56}, with a few modifications. Briefly, 5 mg of tissue were homogenized in 1X PBS on an Omni bead ruptor (Omni International, Kennesaw, GA) prior to adding 20 μ l of internal solution containing U-¹³C labeled analytes (¹³C₃ sodium lactate, ¹³C₄ succinic acid, ¹³C₄ fumaric acid, ¹³C₄ alpha ketoglutaric acid, ¹³C₄ malic acid, ¹³C₄ aspartic acid, ¹³C₅ 2-hydroxyglutaric acid, ¹³C₅ glutamic acid, ¹³C₆ citric acid, ¹³C₂, ¹⁵N glycine, ¹³C₂ sodium pyruvate). For plasma, 50 μ l were used. The proteins were removed by adding 300 μ l of chilled methanol and acetonitrile solution to the sample mixture. After drying the supernatant in the speedvac, the sample was derivatized with ethoxime and then with MtBSTFA + 1% tBDMCS before it was analyzed on an Agilent 5977B GC/MS (Santa Clara, California) under single ion monitoring conditions using electron ionization. Concentrations of lactic acid (*m/z* 261.2), fumaric acid (*m/z* 287.1), succinic acid (*m/z* 289.1), ketoglutaric acid (*m/z* 360.2), malic acid (*m/z* 419.3), aspartic acid (*m/z* 418.2), 2-hydroxyglutaric acid (*m/z* 433.2), cis aconitic acid (*m/z* 459.3), citric acid (*m/z* 591.4), and isocitric acid (*m/z* 591.4), glutamic acid (*m/z* 432.4) were measured against 7-point calibration curves that underwent the same derivatization procedure using ethoxyamine and MtBSTFA. The injection volume was 1 μ L, and the flow rate was 1 mL/min. Table 8 shows the temperature program. Agilent MassHunter was used for data processing.

Table 8. GC temperature gradient program

	Rate ($^{\circ}$ C/min)	Value ($^{\circ}$ C)	Hold Time (min)	Run Time (min)
Initial		120	0.5	0.5
Ramp 1	25	188	0	3.22
Ramp 2	7	275	0	15.649
Ramp 3	35	320	1.066	18