

West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 1 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019
Issued: 07-26-2019 Valid from: 07-26-2019	Validity area: UC Davis Genome Center, Metabolomics Core and Research Laboratories	
Responsible: Ying Yng Choy	Secondary: Bryan Roberts, Arpana Vaniya	
This SOP supersedes: new	Approved: Oliver Fiehn	

1. Instruments

- Thermo Scientific Q Exactive mass spectrometer
- Thermo Scientific Vanquish UHPLC Systems
- Ultrasonicator
- VWR VX-2500 Multi-tube vortexer
- Eppendorf Centrifuge 5424


2. Chemicals and consumables


- Kinetex PFP column (1.7 μ m, 2.1 mm X 30 mm) (00A-4476-AN)
- Security Guard ULTRA Holder, for UHPLC Columns 2.1 to 4.6mm ID (AJ0-9000)
- SecurityGuard™ ULTRA Cartridges UHPLC PFP 2.1mm ID Columns, 3/Pk (AJ0-8787)
- Pipettes
- Pierce™ LTQ ESI Positive Ion Calibration (88322)
- Pierce™ Negative Ion Calibration Solution (88324)
- Acetonitrile: Optima LC/MS Grade Fisher Chemical (A955-4, 4L)
- Water: Optima LC/MS grade Fisher Chemical (W6-4, 4 L)
- Formic Acid: Fluka Mass Spec Grade (F0507, 250 mL)

3. Procedure

3.1 Pre-run procedures


3.1.1 Instrument calibration


- a. Click on Tune icon . Go to **File > Load tune file**. (i.e. Polyphenol_Tune_07122019.mstune)
- b. Select appropriate polarity and click **Apply**.



- c. Use the Thermo Scientific Pierce ESI Negative Ion Calibration Solution for negative mode calibration and Pierce™ LTQ ESI Positive Ion Calibration for positive mode calibration. Load the ready use reference standard solution into a syringe and inject directly into the ion source using syringe pump and designated peek tubing.

West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 2 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

d. Click on syringe pump settings . Enter flow rate ($\mu\text{L}/\text{minutes}$) (i.e. 5, 10, etc), then click **Apply** and **Start** to begin the calibration procedure.




Syringe type	Hamilton
Volume (μL)	250.0
Syringe inner diameter (mm)	2.303
Flow rate ($\mu\text{L} / \text{min}$)	10.000

Buttons: Apply, Help


Manual control buttons: Manual control, Start, Stop

e. In (-) mode, check the profile of the calibrant and ions of m/z 265.148 (SDS), m/z 514.284 (sodium taurocholate), and Ultramark 1621. The TIC variation should be less than 15%. If the TIC variation is higher than 15%, clean the ion source and repeat the instrument calibration.




West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 3 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

f. Once the signal intensity has stabilized, check **Mass calibration (neg)** and click **Calibrate**. The instruments should be re-calibrated weekly.



g. Repeat steps b-d for positive mode calibration.




h. In (+) mode, check the profile of the calibrant and ions of m/z 195.088 (Caffeine), m/z 524.265 (MRFA) and Ultramark 1621.



i. If the TIC variation is higher than 15%, clean the ion source and repeat the instrument calibration.

West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 4 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

3.1.1 Column equilibration

- Click on Thermo Xcalibur icon .
- From the menu bar, click  for main menu or  for sequence setup.
- To take control of the UHPLC module, click **Status** and then **Direct Control** (bottom).
- This will direct you to the pump module page.
- Purge the system by clicking **More Options**. Make sure appropriate eluent is selected (i.e. A1, B1 or A2, B2). Click 100% on line A and purge for 3 minutes at 3 mL/minutes with 50:50 = (acetonitrile: water). Repeat for line B.
- Then purge line A with mobile phase (A) and line B with mobile phase (B) for 3 minutes at 3 mL/minutes to remove acetonitrile and water.
- Place the column with guard column into the column compartment. Equilibrate the column with mobile phase (A) and mobile phase (B) at ratio of 93:7 (starting condition) by setting the pump flow rate at 0.1 mL/minutes; increase the flow rate to 0.4 mL/minutes gradually over 5 minutes. Leave at this flow rate for 30 minutes.

3.2 Preparation of solutions

- Mobile phase A** (100% LC/MS grade water + 0.1% Formic Acid)
 - Pre-rinse three times 1 L glass bottle with pure acetonitrile.
 - Add exactly 1000 mL LC/MS grade water.
 - Add 1 mL formic acid to the a 1 L glass bottle.
 - Sonicate for 10 minutes at room temperature.
- Mobile phase B** (100% ACN + 0.1% Formic Acid)
 - Pre-rinse three times 1 L glass bottle with pure acetonitrile.
 - Add exactly 1000 mL LC/MS acetonitrile.
 - Add 1 mL formic acid to the a 1 L glass bottle.
 - Sonicate for 10 minutes at room temperature.
- Resuspension solvent
 - Prepare 100 mL of 7% ACN in water in a graduate flask.
 - Aliquot 100 uL of hippuric acid-d5 iSTD (100 ppm), 100 uL of quercetin-d3 iSTD (1000 ppm), and 200 uL for all other iSTDs (100 ppm) (See Table 1 & 2) into a 20 mL volumetric flask.
 - Bring to the final volume of the solution to 20 mL with 7% ACN in water.
 - Sonicate for 10 minutes.
 - Store at 4°C.
- Preparation of samples for analysis
 - Sonicate resuspension solvent for 5 minutes and swirl gently before use.
 - Resuspend samples in 100 µL of the resuspension solvent depending on sample concentration.
 - Vortex sample for 30 seconds.
 - Sonicate samples for 5 minutes.
 - Centrifuge for 2 minutes at 16,100 rcf.
 - Transfer supernatant to LC-MS micro insert in an amber vial.
 - Cap and load into auto sampler.

West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 5 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

3.3 Run sequence setup




- a. From the menu bar, click  for sequence setup.
- b. Create run list on excel spreadsheet and save as .csv.
- c. From the menu bar, chose **File > Import sequence > .csv file**
- d. After running the pre-run sequence (3 x blank injection), inject the following:
 - i) Method Blank
 - ii) BioRec EDTA plasma/ NIST SRM3291
 - iii) Pool samples
 - iv) 10 samples
 - v) Repeat steps i-iv

Table 1 Internal standards for positive mode

<u>Name</u>	<u>Concentration in Resuspension Solvent (ppm or ug/mL)</u>	<u>MS1 (m/z)</u>	<u>RT (min)</u>
Hippuric acid-d5 iSTD	0.5	185.097	0.76
Caffeine-d10 iSTD	1	204.144	0.77
(-)-Epigallocatechin Gallate-d3/d4 iSTD	1	462.111	1.44
trans-Cinnamic Acid -d5 iSTD	1	154.091	3.01
trans-Resveratrol-d4 iSTD	1	233.111	3.03
Daidzein-d4 iSTD	1	259.09	3.16
Quercetin-d3 iSTD	5	306.069	3.92
Genistein-d4 iSTD	1	275.085	4.32
Apigenin-d5 iSTD	1	276.092	4.62
2-hydroxyfluorene-d9 iSTD	1	192.136	5.39
CUDA iSTD	1	341.28	5.71
Reserpine-d9 iSTD	1	618.337	7.58

West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 6 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polyphenols07262019

Table 2 Internal standards for negative mode

<u>Name</u>	<u>Concentration in Resuspension Solvent (ppm or ug/mL)</u>	<u>MS1 (m/z)</u>	<u>RT (min)</u>
Hippuric acid-d5 iSTD	0.5	183.081	0.76
(-)-Epigallocatechin Gallate-d3/d4 iSTD	1	460.095	1.37
trans-Cinnamic Acid -d5 iSTD	1	152.075	2.95
trans-Resveratrol-d4 iSTD	1	231.095	2.97
Daidzein-d4 iSTD	1	257.075	3.11
Quercetin-d3 iSTD	5	304.053	3.86
Genistein-d4 iSTD	1	273.07	4.27
Apigenin-d5 iSTD	1	274.076	4.59
2-hydroxyfluorene-d9 iSTD	1	190.121	5.38
CUDA iSTD	1	339.265	5.71
Betulinic acid-d3 iSTD	1	458.372	6.18
Reserpine-d9 iSTD	1	616.322	7.54

3.4 Reverse Phase analysis method

a. The methods for reverse phase analysis are located under the folder C:\Xcalibur\methods\2019\CORE\Polyphenol

Positive ion mode: RP_pos_CE203040_A2B2_10min


Negative ion mode: RP_neg_CE203040_A2B2_10min

Polarity switching mode: RP_pos_neg_CE203040_A2B2_10min

b. The acquisition method for the reverse phase analysis method is as shown below:


3.4.1 Positive ion mode

i) Pump




West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 8 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

ii) Sampler Module




iii) Column Compartment




West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 9 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

iv) System




v) Startup Shutdown



West Coast Metabolomics Center	SOP Standard Operating Procedure	Page 11 of 10
Date: 07-26-2019	Polyphenol Analysis by UPLC-HRMS-Orbitrap	Polypenols07262019

3.4.2 Negative ion mode

The parameters that vary from the positive mode are the following:



3.5 Column storage

Use this procedure to avoid precipitation mobile phase additives on the column and in the system.

- Flush column with 50:50 (acetonitrile:water) by setting the pump flow rate to 0.1 mL/minutes and increase the flow rate to 0.4 mL/minutes over 5 minutes; keep the column at this flow rate for 10 minutes.
- Flush column with 100% acetonitrile by setting the pump flow rate to 0.1 mL/minutes and increase the flow rate to 0.4 mL/minutes over 5 minutes; keep the column at this flow rate for 10 minutes.
- Remove the column from the system.
- Store the column in the box until the next batch analysis.