

Sample preparation

Extraction solvent: MAA(1:1:1) - Methanol : Acetonitrile : Acetone with internal standards (1:25 internal standard mixture in extraction solvent).

Reconstitution solvent: Methanol:H₂O (2:98)

Procedure:

- Pipette 50 μ L of serum from the sample vials into a labeled 1.5 ml micro-centrifuge tube
- Add 200 μ L of extraction solvent to all tubes.
- Vortex for 5 min, then let sit 30 min at 4°C. Vortex again and let sit at -20°C for 1 h, then centrifuge 10 min at 15,000 rpm.
- Transfer 200 μ L of supernatant to labeled labelled 1.5 ml micro-centrifuge tube.
- Dry all samples under nitrogen stream at room temperature.
- Reconstitute in 80 μ L of reconstitution solvent.
- Vortex for 5 min, centrifuge 5 min at 15,000 rpm.
- Transfer supernatant to autosampler vials with inserts, discard pellets.