

Metabolite extraction for HepG2 cells

1. Culture medium of HepG2 cells in 100-mm dish was aspirated and cells were washed twice by 5% mannitol solution (10 mL first and then 2 mL).
2. The cells were then treated with 800 μ L of methanol and left at rest for 30 s in order to inactivate enzymes.
3. The cell extract was then treated with 550 μ L of Milli-Q water containing internal standards (H3304-1002, Human Metabolome Technologies (HMT), Tsuruoka, Yamagata, Japan) and left at rest for another 30s.
4. The extract was obtained and centrifuged at $2,300 \times g$ and 4°C for 5 min and then 800 μ L of upper aqueous layer was centrifugally filtered through a Millipore 5-kDa cutoff filter (UltrafreeMC-PLHCC, HMT) to remove macromolecules ($9,100 \times g$, 4°C , 120 min).
5. The filtrate was centrifugally concentrated and re-suspended in 50 μ L of Milli-Q water for metabolome analysis at HMT.