


UR47
PURGE
300K



Current Data Parameters
NAME urine-paper-160923
EXPNO 30
PROCNO 1

F2 - Acquisition Parameters
Date_ 20160923
Time 14.39 h
INSTRUM spect
PROBHD Z117768_0063 (
PULPROG zgpurge-gprefocus-4
TD 65536
SOLVENT urine
NS 32
DS 4
SWH 12019.230 Hz
FIDRES 0.366798 Hz
AQ 2.7262976 sec
RG 87.01
DW 41.600 usec
DE 42.44 usec
TE 300.0 K
D1 4.0000000 sec
D16 0.0002000 sec
D20 0.0002000 sec
D21 0.0002000 sec
TD0 1
SFO1 600.1328246 MHz
NUC1 1H
P1 20.07 usec
P2 40.14 usec
PLW1 5.55840015 W
PLW9 0.00008953 W
GPNAM[1] SMSQ10.100
GPZ1 -13.17 %
GPNAM[2] SMSQ10.100
GPZ2 52.68 %
GPNAM[3] SMSQ10.100
GPZ3 -17.13 %
GPNAM[4] SMSQ10.100
GPZ4 68.52 %
P16 1000.00 usec



F2 - Processing parameters
SI 65536
SF 600.1299592 MHz
WDW EM
SSB 0
LB 0.30 Hz
GB 0
PC 0